

ICHTHOPASTE[◇] VISCOPASTE[◇] PB7 ZIPZOC[◇]

Effective lower limb care for your patients

The management of venous leg ulceration and the related chronic skin conditions can be a stressful and painful experience for the patient and time consuming for the clinician. Where venous insufficiency exists, compression therapy remains the cornerstone of treatment, however this may be combined with paste bandages to aid treatment, enabling cost effective care for up to 7 days and reducing the number of dressing changes¹.

Viscopaste
10% Zinc Oxide



ZipZoc
20% Zinc Oxide

IchtHopaste
6.32 % Zinc Oxide
2% Ichthammol

NHS Catalogue
IchtHopaste 7.5cm x 6m
Viscopaste 7.5cm x 6m
ZipZoc Stocking

FP10 Drug Tariff
IchtHopaste 7.5cm x 6m
Viscopaste 7.5cm x 6m
ZipZoc Stocking (4's)
ZipZoc Stocking (10's)

Order Code
EFA 051
EFA 011
EFA 029

PIP Code
033-2668
033-2734
234-3739
234-3747

Evolan Pharma AB
For product information please visit www.evolan.se or contact us on
e-mail, pastebandages@evolan.se tel.+44 (0)7554133321

EVOLAN.

1. Williams C. Examining the range of medicated and paste-impregnated bandages. British Journal Of Nursing 1999; 8(15): 119-120

ZIPZOC (zinc impregnated medicated stocking) **Abbreviated prescribing information.** Refer to summary of product characteristics before prescribing.
Composition: A sterile rayon stocking impregnated with an ointment containing 20% Zinc Oxide. Each stocking (80cm x 7cm) contains about 41.5g of ointment. **Therapeutic indications:** Aid to the treatment of chronic leg ulcers and other associated skin conditions. Where chronic venous insufficiency exists, the medicated stocking can be used as a primary contact layer under compression bandaging or hosiery. **Method of Administration:** Because ZIPZOC is preservative free, it must be changed no less frequently than weekly. For topical application, ZIPZOC should be applied to cover the lower leg from the base of the toes to below the knee. All folds should be smoothed out. To protect clothing, a suitable outer bandage should be worn. If chronic venous insufficiency exists, ZIPZOC may be used as a primary contact layer under compression therapy and may be left in situ for one week, before re-application.
Contra-indications: Arterial leg ulcers. Known hypersensitivity to zinc oxide or the excipients in the ointment. **Special Warnings and Special Precautions:** None known. **Interactions:** None known. **Pregnancy and Lactation:** There is no evidence to suggest ZIPZOC should not be used during pregnancy. **Undesirable Effects:** The skin of leg ulcer patients is easily sensitised to some topical medications. ZIPZOC contains no preservatives thereby reducing the risk of skin reactions. Reported reactions include, rash, erythema, itching and maceration of the wound edge. If the outer bandage is inappropriate or the medicated stocking is not changed with sufficient frequency it can cause the stocking to dry out. **Excipients:** Liquid Paraffin Ph. Eur. White Soft Paraffin BP. **Incompatibilities:** None reported. **Shelf Life:** 3 years. **Storage:** At or below 30°C. Keep all medicines out of the reach of children. **Packaging:** Polyethylene aluminium foil laminated pouches, each pouch containing a single medicated stocking. The outer carton will contain either 4 or 10 pouches. **Legal Category:** P **Market authorization number:** PL44616/0002 **Date of first authorisation/renewal of the authorisation:** 01 October 2006 **Date of revision of text:** 26 October 2018 **Trade Mark of Evolan Pharma AB** **Date of preparation:** October 2018